cálido confort

Termómetro infrarrojo
infrared thermometer
Termômetro infrared
Thermomètre infrarouge

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ENGLISH

WARNING

Please read the instructions for use carefully prior to using the product. Store these in a safe place for future reference.

Children should always use this product under supervision, this product is not a toy.

This is solely intended for domestic use.

In order to prevent any risk of danger, do not handle the product. Only qualified technical personnel from the brand's Official Technical support service may carry out repairs or procedures on the device.

Any incorrect use or improper handling of the product shall render the warranty null and void. Only an authorised technical support service centre may carry out repairs on this product.

B&B TRENDS, SL. disclaims all liability for damages that may occur to people, animals or objects, for the non-observance of these warnings.

PARTS AND COMPONENTS:

- LCD screen
- Temperature sensor lens and tip
- Tip for the forehead thermometer mode
- Disposable ear sheath
- Protection cover
- On/memory button
- Temperature scan button
- Battery compartment *

* Batteries not included.

INSTRUCTIONS FOR USE

Place 2 AAA-size batteries in the appliance. Ensure that they are placed in accordance with the polarity indicated in the appliance. Close the battery compartment.

The thermometer should be left in the location where the reading is to take place at least 30 minutes prior to this in order to avoid

temperature changes in the product.

A disposable sheath (D) is necessary for the middle ear reading mode, and the tip for the forehead thermometer mode (C). Otherwise, the apparatus may not take the temperature reading and the symbol will flash on the LCD.

In order to avoid deviations in the temperature reading, ensure that the ear and forehead are clean.

Any scratches on the membrane of the disposable sheath may lead to inaccurate readings.

If the LCD display does not change when you press any of the buttons, remove the batteries from the battery compartment and then insert them again.

Do not use the thermometer to take the ear temperature reading if you have any inflammation or injury to this area.

We recommend taking a reading on a normal day when you feel well in order to have a reference of the regular temperature. This allows you to see the difference.

Keep in mind that the forehead and ear readings do not necessarily have to display the same result, since there can be variations between one reading and another.

Steps for taking a reading:

- Press the on/off button (E).
- The LCD screen will turn blue displaying all the colour values (Fig. 1) and you will hear an alarm. After this initial check, the display will show the temperature 0.0 indicating that it is now ready to take a reading.

Select reading unit (°C/ °F)

1. Press and hold down the on/off button until the display shows C-F° (Fig. 2). While this image is displayed, release the on/off button to change from Celsius to Fahrenheit or vice versa.
2. To change the reading unit again, proceed as above.

Select the reading mode (ear/front):

1. Press and hold down the on/off button. The screen will initially display the option to select the reading unit. If you keep it held down, you will see three lines and the current reading mode (Fig. 3). While this image is displayed, release the on/off button to change from one mode to another.
2. To change the reading unit again, proceed as above.

EAR reading mode:

Ensure that the disposable sheath is firmly fitted onto the product and that the symbol does not appear on the screen .

In order to prevent the transmission of germs, use a new sheath every time you use the product. If the sheath or its membrane tear, discard and replace with a new one.

Straighten up the ear canal by tugging the ear slightly upwards and backwards. This allows both the tip and lens of the temperature sensor to be directed to the eardrum (Fig. 4).

Slowly insert the tip into the ear canal and press the temperature scan button for 1 second to start the reading. When you hear a beep, this means that the reading has been completed.

ENGLISH

Remove the thermometer from the ear.
The LCD display will show the result of the reading in the the reading unit you have selected in advance.

FOREHEAD reading mode:

Ensure that the forehead reading mode tip is firmly fitted onto the product and that the following symbol does not appear on the screen .

Place the reading tip perpendicular to the forehead by resting against it (Fig. 5) Taking the reading at the centre of the forehead. Press the temperature scanning button for 1 second to start the reading. When you hear a beep, this means that the reading has been completed. Remove the thermometer from the forehead.
The LCD display will show the result of the reading in the the reading unit you have selected in advance.

The instrument will automatically turn off after one minute. If you wish to turn it off manually, press and hold down the on/off button until the display shows the word “OFF” .

Repeat the reading if you do not trust the validity of the reading values. In the event that you believe that the results are not correct, take the appliance to a Technical Support Service.

As this is a home use product, it does not need to be calibrated. The readings of this product are indicative only and you should check any high temperature with your doctor.

Review memory of previous readings:

Press the on/off button. When you press the on/off button, the screen will firstly display the memory (from 10 to 1), followed by the result of this reading.

TROUBLESHOOTING GUIDE:

Possible malfunctions:

PROBLEM		SOLUTION
When pressing the on/off button, nothing appears on the screen.	<ul style="list-style-type: none"> • The batteries have expired. • The polarity with which the batteries were placed is not correct. • The batteries are not making proper contact. 	<ul style="list-style-type: none"> • Change the batteries with new ones. • Insert the batteries properly. • Wipe the cavity where the batteries sit with a dry cloth.
The LCD screen displays the word Err	<ul style="list-style-type: none"> • The thermometer cannot perform readings when the room temperature falls outside of the range 10 °C to 40°C. 	<ul style="list-style-type: none"> • The thermometer should be left in the location where the reading is to take place at least 30 minutes prior to the same.
The batteries expire quickly	<ul style="list-style-type: none"> • The batteries are faulty 	<ul style="list-style-type: none"> • Use alkaline batteries from well-known manufacturers
The following appears on-screen: 	<ul style="list-style-type: none"> • The batteries are about to expire 	<ul style="list-style-type: none"> • Replace the batteries.

PRODUCT CONSERVATION

Remove the batteries when you do not intend to use the product for an extended period of time to prolong its shelf-life.

Protect the unit against moisture, sudden changes of temperature and solar radiation.

Do not use the appliance close to strong magnetic fields.

Clean with a damp cloth and dry afterwards. Clean every 3-6 months.

Do not immerse the product in water or other liquid.

Ensure that the product is always stored with its protective cover correctly positioned.

Use a swab or a soft cloth to clean the temperature sensor lens, that can be wet with alcohol or hot water.

Conduct regular reviews of the different parts that make up the product and take it to an authorised Technical Support Service in the event of any damage.

ENGLISH

TECHNICAL DATA

MEASURING RANGE	SOLUCIÓN
Reading method	Ear Mode: 34.0°C~44.0°C (93.2°F~111.2°F) Forehead mode: 32.0°C~42.9°C (89.6°F~109.2°F) In the event that the reading temperature is lower than that of the measurement range, the screen will show the word "Lo " In the event that the reading temperature is higher than the range of measurement, the display will show the word " Hi "
Precision of the indication	Ear Mode: +/- 0.2°C (35.0°C~42.0°C) +/-0.4°F (95.0°F~107.6°F) +/-0.3°C (<35.0°C/>42.0°C) +/-0.5°F (<95.0°F/>107.6°F) Forehead mode: +/-0.3°C (32.0°C~42.9°C) +/-0.5°F (89.6°F~109.2°F)
Memories	10
Operating Conditions	From +10 °C to 40 °C Humidity less than 85%
Storage Conditions	From -25 °C to 55 °C Humidity less than 93%
Power	2 X 1.5V batteries (not included)

PRODUCT DISPOSAL

Do not dispose of the product at the end of its shelf-life. Contact the organic collection centre closest to your home in accordance with the waste electrical and electronic equipment (WEEE) Directive 2002/96/EC to contribute to the protection of the environment. Request detailed information regarding their current methods of disposing of packaging and scrappage of used devices from your distributor.

We hope that you will be satisfied with this product.

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Carefully read the instructions for use.

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The EC marking certifies that complies with directive 93/42/EEC.

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Fabricant.

Representante autorizado en la UE.

Authorized Representative in the EU.

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Elimination of waste according to the European directive 2002/96/EC on waste electrical and electronic equipment.

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This appliance complies with directive 93/42/EC on healthcare products.

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Hangzhou Hua'an Medical & Health Instruments Co., Ltd.
Building 2, 1# Fuzhu Nan RD, Wuchang Town, Yuhang
District, Hangzhou, Zhejiang 310023 China

EU representative: Shanghai International Holding Corp.
GmbH (Europe)
EiffestraBe 80, 20537 Hamburg, GERMANY

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C/Catalunya, 24 P.L Ca N'Oller - 08130 Santa Perpetua de Mogoda
Barcelona Spain

WARRANTY REPORT

B&B TRENDS, SL. guarantees compliance of this product for the purpose for which it is intended for a period of two years the entire Spanish territory.

In the case of breakdown during the term of this warranty, you are entitled to repair or to the replacement of the product at no charge if the former is unavailable, unless one of these options prove impossible or disproportionate. In this case, you can then choose a reduction in price or cancellation of the sale, which must be dealt with directly with the sales vendor. This also covers replacement of spare parts provided that the product has been used according to the recommendations specified in this manual for both cases, and has not been tampered with by any third party that is unaccredited by B & B TRENDS, SL.

This warranty does not affect your rights as a consumer in accordance with the provisions of Law 23/2003 dated 10 July 2003 on warranties on the sale of consumer goods in the Spanish territory.

USE OF WARRANTY

Consumers should contact a technical service centre authorized by B & B TRENDS for repair of the product as any tampering thereof by persons not authorized by Daga, as well as misuse or inappropriate use, will lead to the cancellation of this warranty.

The warranty must be fully completed and delivered along with the receipt or delivery docket for the effective exercise of rights under this warranty.

Claims **must take place within 2 years** from deliver of the project according to the provisions of the warranty, and the consumer must inform the sales vendor of the lack of conformity, provided that such claim is made within the statutory period of warranty. This warranty, as well as the invoice, receipt, or the delivery order must be retained by the user to facilitate the exercise of these rights.

For technical service and after-sales care outside the Spanish territory, please submit your query to the point of sale where you purchased the item or by contacting **B&B TRENDS, SL.**

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