cálido confort

## ET-250

Electroestimulador digital  
Digital electrostimulator  
Eletroestimulador digital  
Eletroestímulateur numérique

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## APPLIANCE DESCRIPTION

1. Channel intensity indicator 1
2. Screen
3. Channel intensity indicator 2
4. Keys for treatment for specific areas.
5. Keys for frequency configuration.
6. Keys for choosing the treatment mode
7. Channel intensity regulator 2
8. Channel intensity regulator 1
9. Cover for the battery compartment.
10. Indicator of the treatment zone.
11. A specific area has been chosen.
12. A specific type of treatment has been chosen.
13. Treatment type indicator: Conventional massage
14. Treatment type indicator: Tapping
15. Treatment type indicator: Shiatsu
16. Timer.
17. Selected frequency indicator.

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## IMPORTANT WARNING

**Carefully read the instructions before using the product. Store these in a safe place for future reference.**

## SAFETY MEASURES

### Repair:

Only specialised personnel should perform repairs to this electrical appliance. Any incorrectly performed repair can lead considerable hazards for the user. Always contact the technical assistance service authorised by B&B Trends, S.L for repair, where failure to observe these instructions will lead to the total loss of the product warranty.

### Contraindications:

In no way should this device be used as a substitute for medical treatment.

Children should be supervised by an adult to avoid using the device as a toy.

This device is not designed to be used by persons (including children) with reduced physical, sensory or psychological capacities, or who lack the necessary knowledge, unless they are supervised by a person responsible for their safety.

### Avoid using the device in the following cases:

- On the abdominal area during pregnancy.
- While you are asleep.
- On animals.
- While driving or using heavy machinery.
- After ingesting substances that alter your perception (e.g. analgesic medicine, alcohol, drugs).
- While carrying out activities that could have an unpredictable reaction, becoming dangerous.
- High fever.
- If you suffer from thrombosis or thrombophlebitis

When in any doubt over whether the electrical stimulation device is appropriate for you, consult your doctor. **Consult your doctor before using the device if:**

- You have a pacemaker or electronic implants.
- You suffer from the following diseases or problems: circulatory problems, diabetes, thrombosis, varicose veins, open wounds, contusions, tears on the skin, vein inflammation.
- You suffer from hypertension.
- You suffer from heart disease.
- You suffer from illnesses related to hypersensitivity in the skin.
- You suffer from epileptic fits.

## ENGLISH

### Use:

- Only use the device for the purposes indicated in this instruction manual.
- Do not leave the appliance unattended when it is in operation.
- Do not place the electrodes in areas with erythema, dermatitis or varicose veins.
- Do not place the electrodes in your mouth, nasal cavities or other internal parts of the body.
- Do not place the electrodes on the head or neck areas.
- Do not place the electrodes on the genital area.
- **IMPORTANT:** Do not place the electrodes on areas close to the heart.
- Do not place the electrodes on metallic objects. Take special care with objects such as belts, jewelry or piercings.
- Always use the device with the cables and electrodes provided. Do not introduce any other element to the device.
- Only DAGA electrodes can be used.
- When you apply the treatment, you should make sure at all times to remain distant from short wave and microwaves that can cause interferences.
- It is not permitted to use this device while connecting yourself to high frequency surgical equipment.
- Stop using the device if you feel nausea, dizziness, fever, muscular pain or any other type of discomfort. Consult your doctor if this happens.
- Do not use the device in water or around excessively humid surroundings such as saunas, showers or Jacuzzis.
- The electrodes must be placed on clean and totally dry skin.
- Do not pull the cables to disconnect the electrodes.
- Do not remove the electrodes from the skin while the device is in operation. If you wish to change their position or simply remove them, turn off the device prior to be able to carry out this action.
- Do not make any modifications on the equipment or on any of its parts (including electrodes and cables).
- Do not apply gels or creams on the areas of the skin onto which the electrodes are placed. The electrodes are made from a specific material for applying electrical impulses homogenously onto the skin, without producing current spikes on specific areas. The use of gels or creams during the treatment can alter these characteristics and cause unpleasant sensations for The user.

**B&B TRENDS SL** accepts no liability for any damage that may be caused to people, animals or objects owing to the failure to observe this warning correctly.

### GENERAL INFORMATION:

First check that the appliance is complete and undamaged. Contact the store where you purchased the product in the event of any doubt.

Through your purchase you have acquired:

- 1 Daga ET-250 Electro stimulator.
- 2 Connector cables.
- 4 electrodes (4x4cm).
- 1 Instruction manual

14 \* Batteries not included.

## STIMULATION SIGNS

### Configuration of parameters:

At any point in the treatment the intensity of the impulses can be increased or decreased as well as their frequency.

The intensity is related to the tension in each instance in which it is applied to the muscular fibres. This value depends on different factors (e.g. thickness of the skin, distance between electrodes, muscle mass index or predominant type of muscular fibres in the area). Most of all you should consider that the greater muscle mass in the intermediate area between the electrodes requires a greater intensity from the impulses.

As for the frequency, this is related to the speed at which the massages or muscle stimulations seem to be carried out. The excitation intervals are produced more often at a higher frequency.

Both the frequency values and the intensity values should be configured according to their own sensations. A tingling feeling indicates that the intensity is appropriate and on the other hand, a sensation of discomfort or pain indicates that the intensity should be reduced.

### Output tension:

The maximum voltage that this electro-stimulator provides for a charge of  $500\Omega$  is 15V and occurs when intensity level 10 is selected. On the other hand, when level 1 is selected, the output voltage reduces to 0.15V. The minimum value represents 1% of the maximum value for this product.

The variations in the output voltage occur directly proportional to the changes in charge (at a greater charge a greater voltage and a lower charge and a lower voltage). It should be mentioned that, the electrical current (number of electrons which are made to circulate through the body per unit of time) does not seem to be affected by the charge where it is applied (amount of muscular tissue that it goes through).

The greatest voltage spikes are given when level 10 is selected and the charge is excessively large. In these circumstances the tension can reach 200V (open circuit). They are short spikes in voltage with those which transmit an efficient value (RMS) lower than 2.7V.

### Massage rollers:

Using this type of device it is possible to apply three types of very different massages.

**1. Conventional massage:** A group of trains of impulses are different that remain for 4.5 seconds when each one is activated and then deactivate for 2.6 seconds. This process is repeated frequently. Each train of impulses is made made up of signals that are activated between 100 and 240  $\mu$ s with a frequency of 16 ms.

**2. Tapping Treatment:** A group of discreet impulses of about 0.4 ms are generated and are repeated every 250 ms.

**3. Shiatsu Treatment:** A group of trains of impulses are different that remain for 4 seconds when each one is activated and then deactivate for 1 second. This process is repeated frequently. Each train of impulses is made made up of signals that are activated between

## ENGLISH

100 and 240  $\mu$ s with a frequency of 16 ms.

The choice of the type of massage depends on the preferences of the user.

### Specific Treatments:

Through a combination of the basic configurations mentioned on the previous section a total of 5 specific pre-configured treatments are made available to the user for each one of the main areas of the body. In this way you only have to choose the area you wish to work by pressing the appropriate key and enjoy.

The cycle of the work of each on one of the specific treatments is the following:

#### a. Abdomen:

1. Tapping. Duration 10s. Frequency every 1s.
2. Tapping. Duration 10s. Frequency every 500ms.
3. Tapping. Duration 10s. Frequency every 100ms.
4. Shiatsu. Duration 30s.
5. Conventional. Duration 30s.
6. Tapping. Duration 10s. Frequency every 100ms.
7. Conventional. Duration 30s.
8. Shiatsu. Duration 30s.

#### b. Shoulders:

- 
1. Tapping. Duration 10s. Frequency every 1s.
  2. Tapping. Duration 10s. Frequency every 500ms.
  3. Tapping. Duration 10s. Frequency every 250ms.
  4. Tapping. Duration 10s. Frequency every 100ms.
  5. Tapping. Duration 10s. Frequency every 1s.
  6. Tapping. Duration 10s. Frequency every 250ms.
- 

#### c. Joints:

1. Conventional. Duration 30s.
2. Tapping. Duration 30s. Frequency every 100ms.
3. Tapping. Duration 10s. Frequency every 500ms.
4. Conventional. Duration 30s.
5. Tapping. Duration 30s. Frequency every 100ms.
6. Shiatsu. Duration 30s.

#### d. Sole of the foot:

1. Conventional. Duration 30s.
2. Tapping. Duration 30s. Frequency every 100ms.
3. Conventional. Duration 30s.
4. Tapping. Duration 30s. Frequency every 500ms.
5. Tapping. Duration 30s. Frequency every 100ms.
6. Conventional. Duration 30s.

#### e. Extremities:

1. Tapping. Duration 30s. Frequency every 100ms.
2. Tapping. Duration 30s. Frequency every 1s.
3. Shiatsu. Duration 30s.
4. Tapping. Duration 30s. Frequency every 100ms.
5. Conventional. Duration 30s.
6. Tapping. Duration 30s. Frequency every 1s.

## INSTRUCTIONS FOR USE

Before beginning to use the product, make sure that there is no apparent damage to its parts including cables and electrodes). In case there is any damage to the device, please contact our aftersales service or one of the authorised distributors.

Follow the instructions that are presented below in order to start up the ET-250 stimulator:

1. Remove the product from its original packaging taking special care not to damage any of its parts. Do not use excessive force to remove the product from the box.
2. Remove the cover (9) located on the rear part of the device and insert AAA 4 batteries (batteries not included). Make sure to insert them with the correct pole.
3. Cover for the battery compartment again.
4. The device has two output channels, independent of applied intensity, which can be used simultaneously or individually. Connect one or two connection cables (depending on whether you wish to use one channel or both) on the device.
5. Connect the electrodes to the cables of the device and stick them onto the skin on the area of treatment. Try to always do it on clean, dry skin, with not metallic elements such as jewellery, piercings or belts.
6. Turn of the device by turning the wheels (7) and/or (8) clockwise. Initially you should adjust the level of intensity to value 1.
7. Select one of the 5 programmes designed for the treatment on specific areas of the body (abdomen, shoulders, joints, extremities of the soles of the feet) or one of the 3 massage modes (conventional, tapping or shiatsu).
8. Gently turn of the device by turning the dials (7) and/or (8) clockwise in order to increase the intensity. You should not exceed the intensity of the treatment. You should adjust it in such a way so that it does not become unpleasant.
9. You can at any time modify the frequency of the treatment (5).
10. The treatment finishes automatically after 15 minutes. You can end it before if you wish by turning the dials (7) and/or (8) anticlockwise until the indication "OFF" is shown.
11. Once the session has finished make sure the the device is turned off completely and remove the electrodes. If the device is not going to be used for a long time, take the batteries out of the device and keep them in a safe place.

## PLACEMENT OF THE ELECTRODES

Place the electrodes in such a way that they follow the direction of the muscle. In case of a specific application in one area, we recommend consulting a Doctor to assess the optimum location for the electrodes.

Bear in mind that the greater the distance of separation of the electrodes, the greater the muscle area is to work. They will also assume that the necessary intensity will be greater than for a smaller muscle area.

## ENGLISH

At the end of this manual you will find a table of recommendations for positioning the electrodes for different parts of the body.

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### RESOLUTION GUIDE: Possible anomalies in the functioning

INCIDENCE	CAUSE	SOLUTION
The device does not turn on.	The product does not have batteries or they are inserted incorrectly.	Insert 4 AAA batteries with the correct pole.
	The batteries have died.	Replace the batteries with new ones
	The on/off dial remains "OFF".	Turn the dial clockwise to turn on the device.
	Damaged equipment.	Contact the Daga technical service.
The stimulation is low or non-existent.	Low intensity level.	Slowly increase the intensity through the dial of the corresponding channel.
	On the screen there is a "P".	Select a specific type of treatment (6) or a dedicated programme (4).
	Badly connected electrodes.	Make sure that the cable is correctly connected to the appliance and to the electrodes.
	Dirty electrodes.	Clean the electrodes with cold water as shown in this manual.
	Electrodes in poor conditions.	Replace the electrodes with new ones. You can acquire them and the normal point of purchase or contact the DAGA customer service centre.
	Short-circuitied electrodes (in contact with one another).	Position these electrodes leaving a greater distance between them.
	The distance between the electrodes is excessive.	Position these electrodes leaving a lower distance between them.
	The batteries are about to run out.	Replace the batteries with new ones

INCIDENCE	CAUSE	SOLUTION
The skin in the treated areas has reddened.	The time during which the treatment has been applied is excessive.	Stop applying the treatment to the reddened zone until the state of the skin returns to normal. In the event that the reddening persists consult a Doctor.
	You suffer from illness or disorders related to hypersensitivity in the skin.	Consult your Doctor to know if you can continue to use this device.

## PRODUCT CONSERVATION

If you do not plan to use the appliance for an extended period of time, you should store it in its original packaging in a clean, cool and dry area. Always take the batteries out of the device when it is not going to be used for a long time.

The electrodes should be kept with the original laminate storage plastic while they are not being used. In this way they will not deteriorate or age prematurely. When the electrodes do not stick to the skin, they should be replaced by new ones. For this you can acquire them and the normal point of purchase or contact the DAGA customer service centre.

## PRODUCT CLEANING

This product can only be cleaned by hand, using a soft cloth slightly moistened with a mild soap. The use of abrasive cleaning products or brushes, solvents, benzine, thinner or alcohol is completely prohibited. Never immerse the appliance or any of its parts, except the electrodes, in water or excessively moisten it in order to clean it. Do not re-use the appliance until it is completely dry.

The electrodes can be cleaned with cold water to eliminate dirt that is stuck to them. Avoid the use of hot water as it could deform the electrodes. The electrodes should be left to dry in the open air having been washed and will not be able to be used again until they are completely dry.

## PRODUCT DISPOSAL

Do not dispose of the product at the end of its shelf-life. Contact the organic collection centre closest to your home in accordance with the waste electrical and electronic equipment (WEEE) Directive 2012/19/EU to contribute to the protection of the environment. Request detailed information regarding their current methods of disposing of packaging and scrappage of used devices from your distributor.

We hope that you will be satisfied with this product.

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**Carefully read the instructions for use.**

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**This product complies with European Directive 2012/19/EU on electrical and electronic devices, known as WEEE (Waste Electrical and Electronic Equipment), that provides the legal framework applicable in the European Union for the disposal and reuse of waste electronic and electrical devices. Do not dispose of this product in the bin. Go to the electrical and electronic waste collection centre closest to your home.**

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**The EC marking certifies that complies with directive 93/42/EEC.**

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**Fabricant.**

Representante autorizado en la UE.

**Authorized Representative in the EU.**

Representante autorizado na UE.

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**B&B Trends, S.L.**

C/Catalunya, 24 P.I. Ca N'Oller. 08130 Santa Perpètua de Mogoda  
Barcelona – España

**Shenzhen Jingkehui Electronic Co., Ltd**

5F Building 12, Hengmingzhu Industrial Park, Xiangxing Road,  
Shajing, Bacan, Shenzhen. China

**MEDevice Union, Ltd**

26 York Street. London, United Kingdom. W1U6PZ

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## WARRANTY REPORT

**B&B TRENDS, SL.** guarantees compliance of this product for the purpose for which it is intended for a period of two years for the entire Spanish territory. In the case of breakdown during the term of this warranty, you are entitled to repair or to the replacement of the product at no charge if the former is unavailable, unless one of these options prove impossible or disproportionate. In this case, you can then choose a reduction in price or cancellation of the sale, which must be dealt with directly with the sales vendor. This also covers replacement of spare parts provided that the product has been used according to the recommendations specified in this manual for both cases, and has not been tampered with by any third party that is unaccredited by **B&B TRENDS, SL.** This warranty does not affect your rights as a consumer in accordance with the provisions of Law 23/2003 dated 10 July 2003 on warranties on the sale of consumer goods in the Spanish territory.

## USE OF WARRANTY

Consumers should contact a technical service centre authorised by **B&B TRENDS, SL.** for repair of the product as any tampering thereof by persons not authorised by **B&B TRENDS, SL.**, as well as misuse or inappropriate use, will lead to the cancellation of this warranty. The warranty must be fully completed and delivered along with the receipt or delivery docket for the effective exercise of rights under this warranty.

Claims must take place within three years from delivery of the product according to the provisions of the warranty, and the consumer must inform the sales vendor of the lack of conformity, provided that such claim is made within the statutory period of warranty. This warranty should be retained by the user as well as the invoice, receipt or the delivery docket to facilitate the exercise of these rights. For technical service and after-sales care outside the Spanish territory, please submit your query to the point of sale where you purchased the item or by contacting **B&B TRENDS, SL.**

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