Flexy-Heat wellness

RM-500

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Massage roll

Rolo massageador

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calidad

control

seguridad

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IMPORTANT WARNING

Read the instructions carefully before using the product. Keep them in a safe place for future reference.

SAFETY MEASURES

This product can be used by children aged 8 years and persons with reduced physical, sensory or mental abilities or lack of experience and knowledge, under the supervision of an adult or if they have been provided with the appropriate guidelines on the safe use of the device and understand the dangers involved.

Children must not play with the appliance.

Children must not perform cleaning and maintenance tasks without supervision.

Repair:

Only specialized personnel should perform repairs to this electrical appliance. Always contact the after-sales service authorized by B&B Trends for repair. Failure to observe these instructions will lead to the total loss of the product warranty.

Contraindications:

This device is intended solely to massage various parts of the body. This should never be used to substitute medical treatment.

Avoid using the massage appliance in the following cases:

- During pregnancy.
- While you are asleep.

ENGLISH

- On children.
- On animals.
- While you are driving.
- Following the intake of substances that alter your perception (e.g. painkillers, alcohol, drugs).
- While performing activities where you may have an unexpected reaction that could be dangerous.

Check with your doctor in the event of any doubt about whether the massage appliance is suitable for you.

Check with your doctor prior to using the massage appliance if:

- You have a pacemaker, artificial limbs or electronic implants.
- You suffer from any of the following diseases or problems: circulatory problems, diabetes, thrombosis, varicose veins, open wounds, contusions, tears on the skin, vein inflammation.
- If you have severe pain in the cervical region.

Use:

- Only use the device for the purposes indicated in this instruction manual. Using this appliance for other purposes shall render the warranty null and void.
- When you remove the plastic packaging from the product, do not leave it within the reach of children. There is a danger of suffocation.
- The maximum recommended duration of use is 15 minutes. After this time, allow the appliance to cool for at least 10 minutes before using it again.
- Do not place the device near a heat source under any circumstances.
- Do not use the appliance and consult your doctor if you feel discomfort or illness.
- Do not use the device near your eyes or other sensitive areas of the body.
- Remove the batteries if the device is not going to be used for a prolonged period of time.
- Keep the batteries out of the reach of children to prevent them from choking.
- When using the device, avoid applying the treatment over an excessively long period on a specific area of the body.

B&B TRENDS SL accepts no liability for any damage that may be caused to people, animals or objects owing to failure to observe this warning correctly.

GENERAL INFORMATION:

First check that the appliance is complete and undamaged. Contact the store where you purchased the product in the event of any doubt.

To apply a massage turn on the massager by pressing the on/off switch on the side of the product and slide it slowly massaging the

area without exerting too much pressure on the area to be treated. It can also be used placing it on the cervical and lumbar area or legs, supporting the body on the product.

You can also place the body part (leg, back) resting on the massager and slowly roll up and down for a pleasant and uniform massage.

The massage roller also allows you to perform yoga as it is a product that can also be used as support such as yoga cubes.

Replacing the batteries (not included)

- Open the cover by the closure.
- Remove the battery compartment slowly by pulling the rope holding it.
- Remove the cover that covers the compartment where the batteries are to be placed.
- Insert two 1.5V AA batteries according to the correct poles indicated on the inside of the compartment.
- Place the battery cover on the product again.
- Change the batteries with new ones if there is decrease in the intensity of the vibration.
- Always change all the batteries at once.

STORAGE AND DISPOSAL OF THE PRODUCT

PRODUCT STORAGE

If you do not plan to use the appliance for an extended period of time, you should store it in its original packaging in a clean, cool and dry area, with no load placed on top of the product. Always remove the batteries from the device before storing it.

PRODUCT CLEANING

Remove the cover carefully and hand wash, do not use bleach. Only the cover can be washed.

PRODUCT DISPOSAL

Do not dispose of the product in the rubbish at the end of its shelf-life. Contact the organic collection centre nearest your home in accordance with the Waste Electrical and Electronic Equipment (WEEE) Directive 2012/19/EU to contribute to the protection of the environment.

Request detailed information regarding their current methods of disposing of packaging and scrapping of used devices from your distributor.

Remove the batteries before dismantling the product. Under no circumstances should you throw used batteries into the rubbish, only into a special waste container, or leave them in a special batteries container.

We hope that you will be satisfied with this product.

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This product complies with European Directive 2012/19/EU on electrical and electronic devices, known as WEEE (Waste Electrical and Electronic Equipment), which provides the legal framework applicable in the European Union for the disposal and reuse of waste electronic and electrical devices.

Do not dispose of this product in the bin, instead going to the electrical and electronic waste collection centre closest to your home.

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The product may contain batteries. Remove them before disposing of the product and dispose in special containers approved for this purpose.

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WARRANTY REPORT

B&B TRENDS, SL. guarantees compliance of this product for the purpose for which it is intended for a period of two years for the entire Spanish territory. In the case of breakdown during the term of this warranty, you are entitled to repair or to the replacement of the product at no charge if the former is unavailable, unless one of these options prove impossible or disproportionate. In this case, you can then choose a reduction in price or cancellation of the sale, which must be dealt with directly with the sales vendor. This also covers replacement of spare parts provided that the product has been used according to the recommendations specified in this manual for both cases, and has not been tampered with by any third party that is unaccredited by **B&B TRENDS, SL.** This warranty does not affect your rights as a consumer in accordance with the provisions of Law 23/2003 dated 10 July 2003 on warranties on the sale of consumer goods in the Spanish territory.

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www.daga.eu

USE OF WARRANTY

Consumers should contact a technical service centre authorised by B&B TRENDS, SL. for repair of the product as any tampering thereof by persons not authorised by B&B TRENDS, SL., as well as misuse or inappropriate use, will lead to the cancellation of this warranty. The warranty must be fully completed and delivered along with the receipt or delivery docket for the effective exercise of rights under this warranty.

Claims must take place within three years from delivery of the product according to the provisions of the warranty, and the consumer must inform the sales vendor of the lack of conformity, provided that such claim is made within the statutory period of warranty. This warranty should be retained by the user as well as the invoice, receipt or the delivery docket to facilitate the exercise of these rights. For technical service and after-sales care outside the Spanish territory, please submit your query to the point of sale where you purchased the item or by contacting B&B TRENDS, SL.

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