cálido confort

BT-200

Báscula de baño
Bathroom scale
Balança de banho
Pèse-personne

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ENGLISH

IMPORTANT WARNING

Read the instructions carefully before using the product. Keep them in a safe place for future reference.

SAFETY MEASURES

- This device is intended solely for domestic use.
- Remove the batteries from the product if it is not being used for a long time.
- Keep the batteries out of the reach of children and people with disabilities.
- Do not move the device while in use.
- Do not exceed the weighing capacity of the device.
- Do not use the device with wet feet.
- Do not store the device in upright position.
- Do not use the device with cracked or broken glass.
- Improper use of the product may involve danger and revoke the warranty and liability.

B&B TRENDS SL. accepts no liability for any damage that may be caused to people, animals or objects owing to failure to observe this warning correctly.

GENERAL INFORMATION:

We would like to thank you for choosing Daga.

First check that the appliance is complete and undamaged. Contact the store where you purchased the product in the event of any doubt.

With your purchase you have acquired:

- 1 Daga scale

INSTRUCTIONS FOR USE

The safety instructions marked at the beginning of this manual must be followed at all times when using this product.

Follow the instructions below in order to switch on the appliance:

1. Remove the product from its original packaging taking special care not to damage any of its parts. Do not use excessive force to remove the product from the box.
2. Remove the cover of the batteries on the back of the scale.
3. Place the batteries in the right place under the cover following the indicated polarity, that is, placing the positive "+" facing out.
4. Put the cover back on.
5. Place the device on a hard, rigid, flat, stable horizontal surface free of vibration.

6. Prepare the apparatus according to the unit weight, kg/lb selecting the weight units kg/lb/st pressing the button on the back of the scale, marked with the word "UNIT"
7. Gently place your feet on the surface of the scale to start the device and wait for the symbol "0.0" to appear on the screen.
8. Place yourself in the centre of the platform and avoid sudden movements, keep still and without touching other objects, until the scale shows your weight. The user must wait while the weight is being calculated. The scale indicates that it has completed the weighing by flashing the final value. two times.
9. To save energy, the device goes into sleep mode after several seconds of inactivity.
10. To return to normal operation, simply place your foot lightly on the platform and wait for the symbol "0.0" to appear on the screen.
11. Replace the unit in a safe and stable place.

TARE FUNCTION:

This product has a TARE function that will enable you to weight your baby, pet or suitcase.

Press the circular icon placed under the LCD display when the scale is off .

The symbol showing that the TARE function is activated will flash on the display .

The adult sign will remain on the display  and "0.0" will be shown. Place yourself in the centre of the platform and avoid sudden movements, keep still and do not touch other objects, until the scale shows your weight. The user must wait while the weight is being calculated. The scale indicates that it has completed the weighing by flashing the final value. two times.

Step off the scale and the symbol  and "0.00" will be shown.

Place yourself centre of the platform with your baby or pet in arms, or that you want to weigh, and the scale will show the weight.

Contraindications:

This product is for domestic use only and not for professional purposes.

Children must be supervised by an adult to prevent them using the appliance as a toy.

This device is not designed to be used by persons (including children) with reduced physical, sensory or psychological capacities, or who lack the necessary experience and knowledge, unless they are supervised by a person responsible for their safety.

Repair:

Only specialised personnel should perform repairs to this electrical appliance. Any incorrectly performed repair may lead to considerable risk for the user. In the event of repair, please always contact our aftersales service. Failure to observe these instructions will lead to the total loss of the product warranty.

TROUBLESHOOTING:

Possible operating anomalies

INCIDENT	CAUSE	SOLUTION
The device does not work.	Device not switched on or misplaced battery.	Put on the device, and wait until your weight is shown or press on it with your feet until "0.0" is shown on the display. Check that the battery is installed with positive polarity "+" facing outwards.
	The device is damaged.	Contact the Daga technical service.
Error messages "Err"	Overload	It means that it has exceeded the maximum weight

PRODUCT CLEANING

This product can only be cleaned by hand, using a soft damp cloth and mild soap. The use of abrasive cleaning products or brushes, solvents, benzene, thinner or alcohol is completely prohibited. Never immerse the appliance or any of its parts in water or excessively moisten it in order to clean it. Do not re-use the appliance until it is completely dry.

PRODUCT DISPOSAL

Do not dispose of the product in the rubbish at the end of its shelf-life. Contact the organic collection centre nearest your home in accordance with the Waste Electrical and Electronic Equipment (WEEE) Directive 2012/19/EU to contribute to the protection of the environment.

Request detailed information regarding the current methods for disposing of packaging and used devices from your distributor.

We hope that you will be satisfied with this product.

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This product complies with European Directive 2012/19/EU on electrical and electronic devices, known as WEEE (Waste Electrical and Electronic Equipment), which provides the legal framework applicable in the European Union for the disposal and reuse of waste electronic and electrical devices.

Do not dispose of this product in the bin, instead going to the electrical and electronic waste collection centre closest to your home.

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The product may contain batteries. Remove them before disposing of the product and dispose in special containers approved for this purpose.

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This symbol means you should be careful and that the product surface is slippery when wet.

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C/Catalunya, 24 P.L Ca N'Oller - 08130 Santa Perpetua de Mogoda
Barcelona Spain

WARRANTY REPORT

B&B TRENDS, SL. guarantees compliance of this product for the purpose for which it is intended for a period of two years for the entire Spanish territory. In the case of breakdown during the term of this warranty, you are entitled to repair or to the replacement of the product at no charge if the former is unavailable, unless one of these options prove impossible or disproportionate. In this case, you can then choose a reduction in price or cancellation of the sale, which must be dealt with directly with the sales vendor. This also covers replacement of spare parts provided that the product has been used according to the recommendations specified in this manual for both cases, and has not been tampered with by any third party that is unaccredited by **B&B TRENDS, SL.** This warranty does not affect your rights as a consumer in accordance with the provisions of Law 23/2003 dated 10 July 2003 on warranties on the sale of consumer goods in the Spanish territory.

USE OF WARRANTY

Consumers should contact a technical service centre authorised by **B&B TRENDS, SL.** for repair of the product as any tampering thereof by persons not authorised by **B&B TRENDS, SL.**, as well as misuse or inappropriate use, will lead to the cancellation of this warranty. The warranty must be fully completed and delivered along with the receipt or delivery docket for the effective exercise of rights under this warranty.

Claims must take place within three years from delivery of the product according to the provisions of the warranty, and the consumer must inform the sales vendor of the lack of conformity, provided that such claim is made within the statutory period of warranty. This warranty should be retained by the user as well as the invoice, receipt or the delivery docket to facilitate the exercise of these rights. For technical service and after-sales care outside the Spanish territory, please submit your query to the point of sale where you purchased the item or by contacting **B&B TRENDS, SL.**

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